

Cosmetic Dentistry Specialising in Natural Prosthetics

A: Unit 2, Stanley Business Park, Prospect Road, Burntwood, WS7 0AE

E: futuresmiles2@aol.com T: 01543 677631

## **Procedures Within The Laboratory For Disinfection Control**

Impressions, models & dental appliances returned to the surgery are non sterile, it is the responsibility of the surgery under the Health & Safety at Work Act to ensure work is sterilised on arrival at the surgery.

Under the Health & safety at Work Act it is the responsibility of the surgery to ensure that impressions, models & dental appliances are sterile before being dispatched to the laboratory, but for our own safety, policies are as follows: Impressions

- Rinse impressions under running tap water prior to casting.
- At the sight of blood, spray with an alcohol based disinfectant and air dry.
- All employees must wear gloves for their own protection.
- All employees are asked for their own safety to be vaccinated in accordance with the recommended blood count against Hepatitis B.

## Models & Dental Appliances.

• Appliances returned to the laboratory are firstly rinsed in running water and then immersed in a Milton solution prior to being rinsed again in water.

This concludes our procedures within the laboratory. Please do not hesitate to contact us if further information is required